



Dipper®

Urinalysis Dipstick Control

- Simulates patient testing
- Liquid, made with human urine
- 18-month refrigerated shelf life from date of manufacture*
- 3 months (or 20 dips) open vial stability when stored at 2°C–8°C*
- Designed for use with most urinalysis reagent strips



Dipper Urinalysis Dipstick Control lets you monitor the performance of visual and instrument readings of urinalysis dipsticks by immersing the dipstick into the control—the same way you test your patient samples. We include values for all dipstick analytes plus microalbumin and creatinine, and qualitative results for hCG early pregnancy detection test methods. It can also be used for confirmatory tests and refractometry.

Dipper Analytes*

Bilirubin, Blood, Creatinine, Glucose, hCG, Ketones, Leukocytes, Microalbumin, Nitrite, pH, Protein, Specific Gravity, Urobilinogen

Catalog

1440-01
1442-61

Contents

Dipper Urinalysis Dipstick Control – 2 Level Set
Dipper Urinalysis Dipstick Control – Level 2

Size

6 x 15 mL
6 x 15 mL

Order now: (310) 536-0006

Visit us at: quantimetrix.com

*Refer to the package insert of currently available lots for specific analyte, stability claims and instrument compatibility.

