



Dropper[®] Plus

Point-of-Care Urinalysis Dipstick Control

- Liquid, made with human urine
- 30 days open vial stability at room temperature for maximum portability*
- 18-month shelf life & open vial stability from date of manufacture when stored at 2°C–8°C*
- Designed for use with most urinalysis reagent strips

Dropper Plus POC Urinalysis Dipstick Control eliminates storage problems and provides maximum portability by providing one month of room temperature stability—perfect for sites without refrigeration. Easy-to-use dropper bottles make dispensing simple. We include values for all dipstick analytes plus microalbumin and creatinine; and qualitative results for hCG early pregnancy detection test methods. It can also be used for confirmatory tests and refractometry.

Dropper Plus Analytes*

Bilirubin, Blood, Creatinine, Glucose, hCG, Ketones, Leukocytes, Microalbumin, Nitrite, pH, Protein, Specific Gravity, Urobilinogen

Catalog #	Contents	Size
1440-04	Dropper Plus POC Urinalysis Dipstick Control – 2 Level Set	10 x 5 mL
1440-06	Dropper Plus POC Urinalysis Dipstick Control – 2 Level Set	2 x 5 mL

Order now: (310) 536-0006

Visit us at: quantimetrix.com

*Refer to the package insert of currently available lots for specific analyte, stability claims and instrument compatibility.

